

Hygo® Jet Pro Bio



Simple handling for everyday clinical use

Hygo Jet Pro Bio is a multi-fit aspirator tube designed for dental teams who value **simplicity and comfort**. The tube is **made from a partly biobased material** derived from renewable sources, supporting a more sustainable clinical practice. Hygo Jet Pro Bio **offers efficient, flexible workflows through its multi-fit design**, ensuring comfortable handling for the clinician and a gentler experience for the patient.

Fits both 11 mm and 16 mm suction systems

Optimized angle for natural hand positioning

Integrated ventilation holes and rounded edges to support patient comfort



Scan for more
information

ORSING

www.orsing.se

Hygo® Jet Pro Bio



Features


Multi-Fit – One tube, two systems

Fits directly with both 11 mm and 16 mm suction systems. The adapter-free design helps create a cleaner setup, fewer loose components and supports a more streamlined workflow.



Ergonomic Design for Everyday Use

The optimised angle is designed to support a natural hand position during treatment, improve access behind the rubber dam and promote comfortable handling during routine procedures.




Designed to support Patient Comfort

Hygo Jet Pro Bio features integrated ventilation holes designed to help minimise soft tissue suction, while the rounded shape is intended to feel gentle during intraoral use. These features are intended to support patient comfort during treatment.



Autoclavable

Autoclavable at 134°C while maintaining its structural integrity. It fits easily into established sterilisation procedures.



Technical specifications:

Length: 138 mm
Fits: Ø 11 mm & 16 mm
Single use

Prod. No: HJPB10

Content: 10 pcs

ORSING
www.orsing.se